




## DNA GEL ELECTROPHORESIS SYSTEM

*Rabrotech horizontal gel electrophoresis system provides the best quality in Innovative design and manufacture, acquired during our existence. Each model Incorporates the many design and safety features, either from within our in house product development team or laboratories worldwide. As a result, all our models are manufactured and finished to the highest specification, easy to use, and with safety always paramount.*



### QR CODE SYSTEM

QR CODE SYSTEM HELPS THE COUSTOMER TO TRACK :  
AUTHENTICITY  
MANUFACTURING DATE  
WARRENTY  
COUSTOMER CARE LINK

### BUILT-IN LEVELLING FEATURE

OUR SPECIAL BUILT IN LEVELING FEATURE HELPS LEVELING GEL TRAY IN ORDER TO RUN  
THE GEL PROPERLY LEVELING GEL TRAY IS IMPORTANT ASPECT OF  
RUNNING GEL IN ELECTROPHORESIS PROCESS  
(LEVELLING FEATURE NOT AVAILAIBLE IN MINI & MINI PLUS)



### GEL CASTER

EASY TO CAST THE GEL WITH SILICON GASKIT  
LEAK PROOF GEL CASTER



## Maxi Horizontal DNA Gel Electrophoresis Unit RE003A

Cat No RE003A MAXI Horizontal DNA Gel Electrophoresis System

DNA Electrophoresis. Gel electrophoresis separates DNA fragments by size in a solid support medium (an Agarose gel.)

|                     |          |
|---------------------|----------|
| Minimum Sample      | 26 wells |
| Maximum sample      | 416wells |
| Gel Size (1)        | 25x22cm  |
| Gel Size (2)        | 20x15cm  |
| Buffer Volume       | 2.8Ltr.  |
| Comb slots distance | 2.5cm    |
| Voltage             | 150-200V |

| Products Code  | Product                                           | Size      | Qty.  |
|----------------|---------------------------------------------------|-----------|-------|
| RE003-UT       | Maxi Plus horizontal gel running tank 35x45x6.5cm | 3Ltr.     | 1     |
| RE003-RT       | Gel Casting Tray with leveller (Gel running try)  | 25X22 cm  | 1     |
| RE003-CT       | Multi Try Gel caster with silicon mat             | 33x27cm   | 1     |
| RE003-PLUS-W   | 0.3mm thick platinum wire                         | 27cm      | 1pair |
| RE15C          | Power leads with 4mm output connectors            | 1.5 meter | 1pair |
| RE003- C16-1   | 16 Well comb                                      | 1mm       | 4     |
| *RE003-C26-1.5 | 26Well comb                                       | 1.5mm     | 4     |

(\*)=multi-channel compatible comb